

## DATASHEET

# CDSN Rabbit Polyclonal Antibody

CAT. NO. APA15259

### KEY FEATURES

Target	CDSN	Source / Host	Rabbit
Reactivity	Human, Mouse, Rat	Clonality	Polyclonal
Applications	WB	Conjugation	Unconjugated
Form / Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.		Storage at -20°C

### BACKGROUND

Important for the epidermal barrier integrity.

### APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:500 - 1:2000
----	----------------

\*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

### PRODUCT OVERVIEW


Description	Rabbit polyclonal antibody to CDSN
Specificity	Recognizes endogenous levels of CDSN protein
Antibody Type	Primary antibody
Immunogen	KLH-conjugated synthetic peptide of human CDSN. The exact sequence is proprietary.
Purification	The antibody was purified by immunogen affinity chromatography.
Molecular Weight	Predicted: 51 kD; Observed: 52 kD
Form/Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Alternative Names	Corneodesmosin; S protein
Gene Symbol	CDSN
Entrez Gene	1041(Human); 386463(Mouse)
SwissProt	Q15517(Human); Q7TPC1(Mouse)

\*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact [info@arexbio.com](mailto:info@arexbio.com) or your local distributor.

\*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

**DATASHEET****CDSN Rabbit Polyclonal Antibody**

CAT. NO. APA15259

**DATA**

Western blot analysis of CDSN expression in A431 (A), NIH3T3 (B), mouse skin (C), rat skeletal muscle (D) whole cell lysates. (Predicted band size: 51 kDa; Observed band size: 52 kDa)

**STORAGE**

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

**NOTE**

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.