

**DATASHEET**

**CstF-77 Rabbit Polyclonal Antibody**

CAT. NO. APA15388

**KEY FEATURES**

Target	CstF-77	Source / Host	Rabbit
Reactivity	Human, Mouse, Rat	Clonality	Polyclonal
Applications	WB	Conjugation	Unconjugated
Form / Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.		Storage at -20°C

**BACKGROUND**

One of the multiple factors required for polyadenylation and 3'-end cleavage of mammalian pre-mRNAs.

**APPLICATION**

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:500 - 1:2000
----	----------------

\*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

**PRODUCT OVERVIEW**


Description	Rabbit polyclonal antibody to CstF-77
Specificity	Recognizes endogenous levels of CstF-77 protein
Antibody Type	Primary antibody
Immunogen	Recombinant fusion protein of human CstF-77. The exact sequence is proprietary.
Purification	The antibody was purified by immunogen affinity chromatography.
Molecular Weight	Predicted: 4; Observed: 77 kD
Form/Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Alternative Names	Cleavage stimulation factor subunit 3; CF-1 77 kDa subunit; Cleavage stimulation factor 77 kDa subunit; CSTF 77 kDa subunit; CstF-77
Gene Symbol	CSTF3
Entrez Gene	1479(Human); 228410(Mouse)
SwissProt	Q12996(Human); Q99LI7(Mouse)

\*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact [info@arexbio.com](mailto:info@arexbio.com) or your local distributor.

\*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

**DATASHEET****CstF-77 Rabbit Polyclonal Antibody**

CAT. NO. APA15388

**DATA**

Western blot analysis of CstF-77 expression in mouse liver (A), rat liver (B) whole cell lysates. (Predicted band size: 4; 12; 82 kD; Observed band size: 77 kD)

**STORAGE**

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

**NOTE**

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.