

DATASHEET

DAZL Rabbit Polyclonal Antibody

CAT. NO. APA15411

KEY FEATURES

| | | | |
|---------------|---|---------------|---------------------|
| Target | DAZL | Source / Host | Rabbit |
| Reactivity | Human, Mouse | Clonality | Polyclonal |
| Applications | WB, IF/ICC | Conjugation | Unconjugated |
| Form / Buffer | Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide. | | Storage at -20°C |

BACKGROUND

RNA-binding protein, which is essential for gametogenesis in both males and females. Plays a central role during spermatogenesis. Acts by binding to the 3'-UTR of mRNA, specifically recognizing GUU triplets, and thereby regulating the translation of key transcripts .

APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

| | |
|--------|----------------|
| WB | 1:500 - 1:2000 |
| IF/ICC | 1:50 - 1:200 |

*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

PRODUCT OVERVIEW

| | |
|-------------------|--|
| Description | Rabbit polyclonal antibody to DAZL |
| Specificity | Recognizes endogenous levels of DAZL protein |
| Antibody Type | Primary antibody |
| Immunogen | Recombinant fusion protein of human DAZL. The exact sequence is proprietary. |
| Purification | The antibody was purified by immunogen affinity chromatography. |
| Molecular Weight | Predicted: 33; Observed: 33 kD |
| Form/Buffer | Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide. |
| Alternative Names | DAZH; DAZL1; DAZLA; SPGYLA; Deleted in azoospermia-like; DAZ homolog; DAZ-like autosomal; Deleted in azoospermia-like 1; SPGY-like-autosomal |
| Gene Symbol | DAZL |
| Entrez Gene | 1618(Human); 13164(Mouse) |
| SwissProt | Q92904(Human); Q64368(Mouse) |

*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact info@arexbio.com or your local distributor.


*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

DATASHEET


DAZL Rabbit Polyclonal Antibody

CAT. NO. APA15411

DATA



Western blot analysis of DAZL expression in SHSY5Y (A), HL60 (B), mouse brain (C) whole cell lysates. (Predicted band size: 33; 35 kD; Observed band size: 33 kD)



Immunofluorescent analysis of DAZL staining in HeLa cells. Formalin-fixed cells were permeabilized with 0.1% Triton X-100 in TBS for 5-10 minutes and blocked with 3% BSA-PBS for 30 minutes at room temperature. Cells were probed with the primary antibody in 3% BSA-PBS and incubated overnight at 4 °C in a humidified chamber. Cells were washed with PBST and incubated with a Alexa Fluor 594-conjugated secondary antibody (red) in PBS at room temperature in the dark.

STORAGE

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

NOTE

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.