

DATASHEET

DHX57 Rabbit Polyclonal Antibody

CAT. NO. APA15455

KEY FEATURES

Target	DHX57	Source / Host	Rabbit
Reactivity	Human, Rat	Clonality	Polyclonal
Applications	WB	Conjugation	Unconjugated
Form / Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.		Storage at -20°C

BACKGROUND

Probable ATP-binding RNA helicase.

APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:500 - 1:2000
----	----------------

*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

PRODUCT OVERVIEW


Description	Rabbit polyclonal antibody to DHX57
Specificity	Recognizes endogenous levels of DHX57 protein
Antibody Type	Primary antibody
Immunogen	Recombinant fusion protein of human DHX57. The exact sequence is proprietary.
Purification	The antibody was purified by immunogen affinity chromatography.
Molecular Weight	Predicted: 49; Observed: 155; 73 kD
Form/Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Alternative Names	Putative ATP-dependent RNA helicase DHX57; DEAH box protein 57
Gene Symbol	DHX57
Entrez Gene	90957(Human)
SwissProt	Q6P158(Human)

*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact info@arexbio.com or your local distributor.

*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

DATASHEET**DHX57 Rabbit Polyclonal Antibody**

CAT. NO. APA15455

DATA

Western blot analysis of DHX57 expression in U87MG (A), rat eye (B) whole cell lysates.
(Predicted band size: 49; 73; 155 kD; Observed band size: 155; 73 kD)

STORAGE

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

NOTE

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.