

DATASHEET

Dermokine Rabbit Polyclonal Antibody

CAT. NO. APA15467

KEY FEATURES

Target	Dermokine	Source / Host	Rabbit
Reactivity	Human, Rat	Clonality	Polyclonal
Applications	WB	Conjugation	Unconjugated
Form / Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.		Storage at -20°C

BACKGROUND

May act as a soluble regulator of keratinocyte differentiation.

APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:500 - 1:2000
----	----------------

*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

PRODUCT OVERVIEW


Description	Rabbit polyclonal antibody to Dermokine
Specificity	Recognizes endogenous levels of Dermokine protein
Antibody Type	Primary antibody
Immunogen	Recombinant fusion protein of human Dermokine. The exact sequence is proprietary.
Purification	The antibody was purified by immunogen affinity chromatography.
Molecular Weight	Predicted: 9; Observed: 47 kD
Form/Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Alternative Names	Dermokine; Epidermis-specific secreted protein SK30/SK89
Gene Symbol	DMKN
Entrez Gene	93099(Human)
SwissProt	Q6E0U4(Human)

*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact info@arexbio.com or your local distributor.

*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

DATASHEET**Dermokine Rabbit Polyclonal Antibody**

CAT. NO. APA15467

DATA

Western blot analysis of Dermokine expression in rat lung (A) whole cell lysates.
(Predicted band size: 9; 14-20; 35-47 kD; Observed band size: 47 kD)

STORAGE

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

NOTE

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.