

## DATASHEET

# FBL14 Rabbit Polyclonal Antibody

CAT. NO. APA15606

### KEY FEATURES

Target	FBL14	Source / Host	Rabbit
Reactivity	Human, Mouse, Rat	Clonality	Polyclonal
Applications	WB	Conjugation	Unconjugated
Form / Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.		Storage at -20°C

### BACKGROUND

Substrate-recognition component of some SCF (SKP1-CUL1-F-box protein)-type E3 ubiquitin-protein ligase complexes. The SCF(FBXL14) complex acts by mediating ubiquitination and subsequent degradation of SNAI1.

### APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:500 - 1:2000
----	----------------

\*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

### PRODUCT OVERVIEW


Description	Rabbit polyclonal antibody to FBL14
Specificity	Recognizes endogenous levels of FBL14 protein
Antibody Type	Primary antibody
Immunogen	Recombinant fusion protein of human FBL14. The exact sequence is proprietary.
Purification	The antibody was purified by immunogen affinity chromatography.
Molecular Weight	Predicted: 45 kD; Observed: 46 kD
Form/Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Alternative Names	FBL14; F-box/LRR-repeat protein 14; F-box and leucine-rich repeat protein 14
Gene Symbol	FBXL14
Entrez Gene	144699(Human); 101358(Mouse)
SwissProt	Q8N1E6(Human); Q8BID8(Mouse)

\*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact [info@arexbio.com](mailto:info@arexbio.com) or your local distributor.

\*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

**DATASHEET****FBL14 Rabbit Polyclonal Antibody**

CAT. NO. APA15606

**DATA**

Western blot analysis of FBL14 expression in HepG2 (A), mouse brain (B), rat thymus (C) whole cell lysates. (Predicted band size: 45 kD; Observed band size: 46 kD)

**STORAGE**

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

**NOTE**

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.