

**DATASHEET**

**GSTT2 Rabbit Polyclonal Antibody**

CAT. NO. APA15776

**KEY FEATURES**

Target	GSTT2	Source / Host	Rabbit
Reactivity	Human, Mouse, Rat	Clonality	Polyclonal
Applications	WB	Conjugation	Unconjugated
Form / Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.		Storage at -20°C

**BACKGROUND**

Conjugation of reduced glutathione to a wide number of exogenous and endogenous hydrophobic electrophiles . Has a sulfatase activity .

**APPLICATION**

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:500 - 1:2000
----	----------------

\*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

**PRODUCT OVERVIEW**


Description	Rabbit polyclonal antibody to GSTT2
Specificity	Recognizes endogenous levels of GSTT2 protein
Antibody Type	Primary antibody
Immunogen	Recombinant fusion protein of human GSTT2. The exact sequence is proprietary.
Purification	The antibody was purified by immunogen affinity chromatography.
Molecular Weight	Predicted: 27 kD; Observed: 27 kD
Form/Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Alternative Names	Glutathione S-transferase theta-2; GST class-theta-2
Gene Symbol	GSTT2
Entrez Gene	2953(Human)
SwissProt	POCG29(Human)

\*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact [info@arexbio.com](mailto:info@arexbio.com) or your local distributor.

\*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

**DATASHEET****GSTT2 Rabbit Polyclonal Antibody**

CAT. NO. APA15776

**DATA**

Western blot analysis of GSTT2 expression in LO2 (A), mouse heart (B), rat liver (C), rat kidney (D) whole cell lysates. (Predicted band size: 27 kD; Observed band size: 27 kD)

**STORAGE**

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

**NOTE**

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.