

DATASHEET

MED25 Rabbit Polyclonal Antibody

CAT. NO. APA16126

KEY FEATURES

| | | | |
|---------------|---|---------------|---------------------|
| Target | MED25 | Source / Host | Rabbit |
| Reactivity | Human, Mouse, Rat | Clonality | Polyclonal |
| Applications | WB, IF/ICC | Conjugation | Unconjugated |
| Form / Buffer | Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide. | | Storage at -20°C |

BACKGROUND

Component of the Mediator complex, a coactivator involved in the regulated transcription of nearly all RNA polymerase II-dependent genes. Mediator functions as a bridge to convey information from gene-specific regulatory proteins to the basal RNA polymerase II transcription machinery. Mediator is recruited to promoters by direct interactions with regulatory proteins and serves as a scaffold for the assembly of a functional preinitiation complex with RNA polymerase II and the general transcription factors. Required for RARA/RXRA-mediated transcription.

APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

| | |
|--------|----------------|
| WB | 1:500 - 1:2000 |
| IF/ICC | 1:50 - 1:200 |

*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

PRODUCT OVERVIEW

| | |
|-------------------|--|
| Description | Rabbit polyclonal antibody to MED25 |
| Specificity | Recognizes endogenous levels of MED25 protein |
| Antibody Type | Primary antibody |
| Immunogen | Recombinant fusion protein of human MED25. The exact sequence is proprietary. |
| Purification | The antibody was purified by immunogen affinity chromatography. |
| Molecular Weight | Predicted: 55; Observed: 78 kD |
| Form/Buffer | Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide. |
| Alternative Names | ACID1; ARC92; PTOV2; Mediator of RNA polymerase II transcription subunit 25; Activator interaction domain-containing protein 1; Activator-recruited cofactor 92 kDa component; ARC92; Mediator complex subunit 25; p78 |
| Gene Symbol | MED25 |
| Entrez Gene | 81857(Human); 75613(Mouse) |
| SwissProt | Q71SY5(Human); Q8VCB2(Mouse) |

*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact info@arexbio.com or your local distributor.


*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

DATASHEET


MED25 Rabbit Polyclonal Antibody

CAT. NO. APA16126

DATA



Western blot analysis of MED25 expression in A549 (A), mouse brain (B), rat kidney (C) whole cell lysates. (Predicted band size: 55; 76; 78; 82; 84 kD; Observed band size: 78 kD)



Immunofluorescent analysis of MED25 staining in U2OS cells. Formalin-fixed cells were permeabilized with 0.1% Triton X-100 in TBS for 5-10 minutes and blocked with 3% BSA-PBS for 30 minutes at room temperature. Cells were probed with the primary antibody in 3% BSA-PBS and incubated overnight at 4 °C in a humidified chamber. Cells were washed with PBST and incubated with a Alexa Fluor 594-conjugated secondary antibody (red) in PBS at room temperature in the dark.

STORAGE

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

NOTE

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.