

DATASHEET

OSTM1 Rabbit Polyclonal Antibody

CAT. NO. APA16384

KEY FEATURES

Target	OSTM1	Source / Host	Rabbit
Reactivity	Human, Mouse, Rat	Clonality	Polyclonal
Applications	WB	Conjugation	Unconjugated
Form / Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.		Storage at -20°C

BACKGROUND

Required for osteoclast and melanocyte maturation and function.

APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:500 - 1:2000
----	----------------

*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

PRODUCT OVERVIEW


Description	Rabbit polyclonal antibody to OSTM1
Specificity	Recognizes endogenous levels of OSTM1 protein
Antibody Type	Primary antibody
Immunogen	Recombinant fusion protein of human OSTM1. The exact sequence is proprietary.
Purification	The antibody was purified by immunogen affinity chromatography.
Molecular Weight	Predicted: 37 kD; Observed: 37 kD
Form/Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Alternative Names	GL; Osteopetrosis-associated transmembrane protein 1; Chloride channel 7 beta subunit
Gene Symbol	OSTM1
Entrez Gene	28962(Human); 14628(Mouse)
SwissProt	Q86WC4(Human); Q8BGT0(Mouse)

*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact info@arexbio.com or your local distributor.

*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

DATASHEET**OSTM1 Rabbit Polyclonal Antibody**

CAT. NO. APA16384

DATA

Western blot analysis of OSTM1 expression in HEK293T (A), A549 (B), mouse brain (C), mouse kidney (D) whole cell lysates. (Predicted band size: 37 kD; Observed band size: 37 kD)

STORAGE

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

NOTE

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.