

DATASHEET

PRPSAP1 Rabbit Polyclonal Antibody

CAT. NO. APA16600

KEY FEATURES

Target	PRPSAP1	Source / Host	Rabbit
Reactivity	Human, Mouse, Rat	Clonality	Polyclonal
Applications	WB	Conjugation	Unconjugated
Form / Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.		Storage at -20°C

BACKGROUND

Seems to play a negative regulatory role in 5-phosphoribose 1-diphosphate synthesis.

APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:500 - 1:2000
----	----------------

*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

PRODUCT OVERVIEW


Description	Rabbit polyclonal antibody to PRPSAP1
Specificity	Recognizes endogenous levels of PRPSAP1 protein
Antibody Type	Primary antibody
Immunogen	Recombinant fusion protein of human PRPSAP1. The exact sequence is proprietary.
Purification	The antibody was purified by immunogen affinity chromatography.
Molecular Weight	Predicted: 39; Observed: 42 kD
Form/Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Alternative Names	Phosphoribosyl pyrophosphate synthase-associated protein 1; PRPP synthase-associated protein 1; 39 kDa phosphoribosypyrophosphate synthase-associated protein; PAP39
Gene Symbol	PRPSAP1
Entrez Gene	5635(Human); 67763(Mouse)
SwissProt	Q14558(Human); Q9D0M1(Mouse); Q63468(Rat)

*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact info@arexbio.com or your local distributor.

*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

DATASHEET**PRPSAP1 Rabbit Polyclonal Antibody**

CAT. NO. APA16600

DATA

Western blot analysis of PRPSAP1 expression in HeLa (A), mouse heart (B) whole cell lysates. (Predicted band size: 39; 42 kD; Observed band size: 42 kD)

STORAGE

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

NOTE

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.