

DATASHEET

Tectonic-2 Rabbit Polyclonal Antibody

CAT. NO. APA17079

KEY FEATURES

Target	Tectonic-2	Source / Host	Rabbit
Reactivity	Human, Mouse, Rat	Clonality	Polyclonal
Applications	WB	Conjugation	Unconjugated
Form / Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.		Storage at-20°C

BACKGROUND

Component of the tectonic-like complex, a complex localized at the transition zone of primary cilia and acting as a barrier that prevents diffusion of transmembrane proteins between the cilia and plasma membranes. Required for hedgehog signaling transduction .

APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:500 - 1:2000
----	----------------

*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

PRODUCT OVERVIEW


Description	Rabbit polyclonal antibody to Tectonic-2
Specificity	Recognizes endogenous levels of Tectonic-2 protein
Antibody Type	Primary antibody
Immunogen	Recombinant fusion protein of human Tectonic-2. The exact sequence is proprietary.
Purification	The antibody was purified by immunogen affinity chromatography.
Molecular Weight	Predicted: 76 kD; Observed: 77 kD
Form/Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Alternative Names	C12orf38; TECT2; Tectonic-2
Gene Symbol	TCTN2
Entrez Gene	79867(Human); 67978(Mouse); 689779(Rat)
SwissProt	Q96GX1(Human); Q2MV57(Mouse); Q3B7D3(Rat)

*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact info@arexbio.com or your local distributor.

*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

DATASHEET**Tectonic-2 Rabbit Polyclonal Antibody**

CAT. NO. APA17079

DATA

Western blot analysis of Tectonic-2 expression in mouse brain (A), mouse lung (B) whole cell lysates. (Predicted band size: 76 kD; Observed band size: 77 kD)

STORAGE

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

NOTE

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.