

DATASHEET

Teneurin-1 Rabbit Polyclonal Antibody

CAT. NO. APA17087

KEY FEATURES

Target	Teneurin-1	Source / Host	Rabbit
Reactivity	Human, Mouse	Clonality	Polyclonal
Applications	WB, IF/ICC	Conjugation	Unconjugated
Form / Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.		Storage at -20°C

BACKGROUND

Involved in neural development, regulating the establishment of proper connectivity within the nervous system. May function as a cellular signal transducer .

APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:500 - 1:2000
IF/ICC	1:50 - 1:200

*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

PRODUCT OVERVIEW

Description	Rabbit polyclonal antibody to Teneurin-1
Specificity	Recognizes endogenous levels of Teneurin-1 protein
Antibody Type	Primary antibody
Immunogen	Recombinant fusion protein of human Teneurin-1. The exact sequence is proprietary.
Purification	The antibody was purified by immunogen affinity chromatography.
Molecular Weight	Predicted: 305 kD; Observed: 280 kD
Form/Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Alternative Names	ODZ1; TNM1; Teneurin-1; Ten-1; Protein Odd Oz/ten-m homolog 1; Tenascin-M1; Ten-m1; Teneurin transmembrane protein 1
Gene Symbol	TENM1
Entrez Gene	10178(Human); 23963(Mouse)
SwissProt	Q9UKZ4(Human); Q9WTS4(Mouse)

*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact info@arexbio.com or your local distributor.


*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

DATASHEET


Teneurin-1 Rabbit Polyclonal Antibody

CAT. NO. APA17087

DATA



Western blot analysis of Teneurin-1 expression in HeLa (A) whole cell lysates. (Predicted band size: 305 kD; Observed band size: 280 kD)



Immunofluorescent analysis of Teneurin-1 staining in C6 cells. Formalin-fixed cells were permeabilized with 0.1% Triton X-100 in TBS for 5-10 minutes and blocked with 3% BSA-PBS for 30 minutes at room temperature. Cells were probed with the primary antibody in 3% BSA-PBS and incubated overnight at 4 °C in a humidified chamber. Cells were washed with PBST and incubated with a AREX® Fluor 594-conjugated secondary antibody (red) in PBS at room temperature in the dark. DAPI was used to stain the cell nuclei (blue).

STORAGE

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

NOTE

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.