

DATASHEET

THSD7A Rabbit Polyclonal Antibody

CAT. NO. APA17104

KEY FEATURES

Target	THSD7A	Source / Host	Rabbit
Reactivity	Human, Mouse	Clonality	Polyclonal
Applications	WB	Conjugation	Unconjugated
Form / Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.		Storage at -20°C

BACKGROUND

Plays a role in actin cytoskeleton rearrangement.

APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:500 - 1:2000
----	----------------

*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

PRODUCT OVERVIEW


Description	Rabbit polyclonal antibody to THSD7A
Specificity	Recognizes endogenous levels of THSD7A protein
Antibody Type	Primary antibody
Immunogen	Recombinant fusion protein of human THSD7A. The exact sequence is proprietary.
Purification	The antibody was purified by immunogen affinity chromatography.
Molecular Weight	Predicted: 185 kD; Observed: 260 kD
Form/Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Alternative Names	KIAA0960; Thrombospondin type-1 domain-containing protein 7A
Gene Symbol	THSD7A
Entrez Gene	221981(Human); 330267(Mouse)
SwissProt	Q9UPZ6(Human); Q69ZU6(Mouse)

*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact info@arexbio.com or your local distributor.

*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

DATASHEET**THSD7A Rabbit Polyclonal Antibody**

CAT. NO. APA17104

DATA

Western blot analysis of THSD7A expression in HEK293T (A) whole cell lysates. (Predicted band size: 185 kD; Observed band size: 260 kD)

STORAGE

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

NOTE

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.