

**DATASHEET**

**TRIM27 Rabbit Polyclonal Antibody**

CAT. NO. APA17177

**KEY FEATURES**

Target	TRIM27	Source / Host	Rabbit
Reactivity	Human, Mouse	Clonality	Polyclonal
Applications	WB, IHC	Conjugation	Unconjugated
Form / Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.		Storage at -20°C

**BACKGROUND**

E3 ubiquitin-protein ligase that mediates ubiquitination of various substrates and thereby plays a role in different processes including proliferation, innate immunity, apoptosis, immune response or autophagy. Ubiquitinates PIK3C2B and inhibits its activity by mediating the formation of 'Lys-48'-linked polyubiquitin chains; the function inhibits CD4 T-cell activation. Acts as a regulator of retrograde transport: together with MAGEL2, mediates the formation of 'Lys-63'-linked polyubiquitin chains at 'Lys-220' of WASHC1, leading to promote endosomal F-actin assembly. Has a transcriptional repressor activity by cooperating with EPC1. Induces apoptosis by activating Jun N-terminal kinase and p38 kinase and also increases caspase-3-like activity independently of mitochondrial events.

**APPLICATION**

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:500 - 1:1000
IHC	1:50 - 1:200

\*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

**PRODUCT OVERVIEW**

Description	Rabbit polyclonal antibody to TRIM27
Specificity	Recognizes endogenous levels of TRIM27 protein
Antibody Type	Primary antibody
Immunogen	Recombinant fusion protein of human TRIM27. The exact sequence is proprietary.
Purification	The antibody was purified by immunogen affinity chromatography.
Molecular Weight	Predicted: 41; Observed: 59 kD
Form/Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Alternative Names	RFP; RNF76; Zinc finger protein RFP; RING finger protein 76; Ret finger protein; Tripartite motif-containing protein 27
Gene Symbol	TRIM27
Entrez Gene	5987(Human); 19720(Mouse)
SwissProt	P14373(Human); Q62158(Mouse)

\*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact [info@arexbio.com](mailto:info@arexbio.com) or your local distributor.


\*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

**DATASHEET**


**TRIM27 Rabbit Polyclonal Antibody**

CAT. NO. APA17177

**DATA**



Western blot analysis of TRIM27 expression in Jurkat (A), mouse testis (B) whole cell lysates. (Predicted band size: 41; 58 kD; Observed band size: 59 kD)



Immunohistochemical analysis of TRIM27 staining in human esophagus formalin fixed paraffin embedded tissue section. The section was pre-treated using heat mediated antigen retrieval with sodium citrate buffer (pH 6.0). The section was then incubated with the antibody at room temperature and detected using an HRP conjugated compact polymer system. DAB was used as the chromogen. The section was then counterstained with haematoxylin and mounted with DPX.

**STORAGE**

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

**NOTE**

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.