

**DATASHEET**

**TRIM47 Rabbit Polyclonal Antibody**

CAT. NO. APA17185

**KEY FEATURES**

Target	TRIM47	Source / Host	Rabbit
Reactivity	Human, Rat	Clonality	Polyclonal
Applications	WB, IF/ICC	Conjugation	Unconjugated
Form / Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.		Storage at -20°C

**BACKGROUND**

E3 ubiquitin-protein ligase that mediates the ubiquitination and proteasomal degradation of CYLD.

**APPLICATION**

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:500 - 1:2000
IF/ICC	1:50 - 1:200

\*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

**PRODUCT OVERVIEW**

Description	Rabbit polyclonal antibody to TRIM47
Specificity	Recognizes endogenous levels of TRIM47 protein
Antibody Type	Primary antibody
Immunogen	Recombinant fusion protein of human TRIM47. The exact sequence is proprietary.
Purification	The antibody was purified by immunogen affinity chromatography.
Molecular Weight	Predicted: 69 kD; Observed: 70 kD
Form/Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Alternative Names	GOA; RNF100; Tripartite motif-containing protein 47; Gene overexpressed in astrocytoma protein; RING finger protein 100
Gene Symbol	TRIM47
Entrez Gene	91107(Human); 217333(Mouse)
SwissProt	Q96LD4(Human); Q8C0E3(Mouse)

\*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact [info@arexbio.com](mailto:info@arexbio.com) or your local distributor.


\*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

**DATASHEET**


**TRIM47 Rabbit Polyclonal Antibody**

CAT. NO. APA17185

**DATA**



Western blot analysis of TRIM47 expression in HeLa (A), rat heart (B) whole cell lysates. (Predicted band size: 69 kD; Observed band size: 70 kD)



Immunofluorescent analysis of TRIM47 staining in L929 cells. Formalin-fixed cells were permeabilized with 0.1% Triton X-100 in TBS for 5-10 minutes and blocked with 3% BSA-PBS for 30 minutes at room temperature. Cells were probed with the primary antibody in 3% BSA-PBS and incubated overnight at 4 °C in a humidified chamber. Cells were washed with PBST and incubated with a AREX® Fluor 594-conjugated secondary antibody (red) in PBS at room temperature in the dark. DAPI was used to stain the cell nuclei (blue).

**STORAGE**

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

**NOTE**

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.