

DATASHEET

Vasorin Rabbit Polyclonal Antibody

CAT. NO. APA17279

KEY FEATURES

Target	Vasorin	Source / Host	Rabbit
Reactivity	Human, Mouse, Rat	Clonality	Polyclonal
Applications	WB	Conjugation	Unconjugated
Form / Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.		Storage at -20°C

BACKGROUND

May act as an inhibitor of TGF-beta signaling.

APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:500 - 1:2000
----	----------------

*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

PRODUCT OVERVIEW


Description	Rabbit polyclonal antibody to Vasorin
Specificity	Recognizes endogenous levels of Vasorin protein
Antibody Type	Primary antibody
Immunogen	Recombinant fusion protein of human Vasorin. The exact sequence is proprietary.
Purification	The antibody was purified by immunogen affinity chromatography.
Molecular Weight	Predicted: 71 kD; Observed: 71 kD
Form/Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Alternative Names	SLITL2; Vasorin; Protein slit-like 2
Gene Symbol	VASN
Entrez Gene	114990(Human); 246154(Mouse)
SwissProt	Q6EMK4(Human); Q9CZT5(Mouse)

*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact info@arexbio.com or your local distributor.

*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

DATASHEET**Vasorin Rabbit Polyclonal Antibody**

CAT. NO. APA17279

DATA

Western blot analysis of Vasorin expression in mouse kidney (A), rat lung (B) whole cell lysates. (Predicted band size: 71 kD; Observed band size: 71 kD)

STORAGE

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

NOTE

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.