

DATASHEET

RGAG4 Rabbit Polyclonal Antibody

CAT. NO. APA17487

KEY FEATURES

Target	RGAG4	Source / Host	Rabbit
Reactivity	Human	Clonality	Polyclonal
Applications	WB, IHC	Conjugation	Unconjugated
Form / Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.		Storage at -20°C

BACKGROUND

retrotransposon Gag like 5. Also known as: retrotransposon Gag-like protein 5; Sushi-Ichi retrotransposon homolog 8; mammalian retrotransposon-derived 5; retrotransposon gag domain containing 4; retrotransposon gag domain-containing protein 4.

APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:500 - 1:1000
IHC	1:50 - 1:200

*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

PRODUCT OVERVIEW

Description	Rabbit polyclonal antibody to RGAG4
Specificity	Recognizes endogenous levels of RGAG4 protein.
Antibody Type	Primary antibody
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the C-terminal region of human RGAG4. The exact sequence is proprietary.
Purification	The antibody was purified by immunogen affinity chromatography.
Molecular Weight	Predicted: 64 kD; Observed: 55 kD
Form/Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Alternative Names	KIAA2001; Retrotransposon gag domain-containing protein 4
Gene Symbol	RGAG4
Entrez Gene	340526(Human)
SwissProt	Q5HYW3(Human)

*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact info@arexbio.com or your local distributor.


*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

DATASHEET


RGAG4 Rabbit Polyclonal Antibody

CAT. NO. APA17487

DATA



Western blot analysis of RGAG4 expression in U251MG (A) whole cell lysates. (Predicted band size: 64 kD; Observed band size: 55 kD)



Immunohistochemical analysis of RGAG4 staining in human brain formalin fixed paraffin embedded tissue section. The section was pre-treated using heat mediated antigen retrieval with sodium citrate buffer (pH 6.0). The section was then incubated with the antibody at room temperature and detected using an HRP conjugated compact polymer system. DAB was used as the chromogen. The section was then counterstained with haematoxylin and mounted with DPX.

STORAGE

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

NOTE

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.