

## DATASHEET

# SALL1 Rabbit Polyclonal Antibody

CAT. NO. APA18239

### KEY FEATURES

Target	SALL1	Source / Host	Rabbit
Reactivity	Human	Clonality	Polyclonal
Applications	WB	Conjugation	Unconjugated
Form / Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.		Storage at-20°C

### BACKGROUND

Transcriptional repressor involved in organogenesis. Plays an essential role in ureteric bud invasion during kidney development.

### APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:500 - 1:1000
----	----------------

\*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

### PRODUCT OVERVIEW


Description	Rabbit polyclonal antibody to SALL1
Specificity	Recognizes endogenous levels of SALL1 protein.
Antibody Type	Primary antibody
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the N-terminal region of human SALL1. The exact sequence is proprietary.
Purification	The antibody was purified by immunogen affinity chromatography.
Molecular Weight	Predicted: 140 kD; Observed: 130 kD
Form/Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Alternative Names	SAL1; ZNF794; Sal-like protein 1; Spalt-like transcription factor 1; Zinc finger protein 794; Zinc finger protein SALL1; Zinc finger protein Spalt-1; HSal1; Sal-1
Gene Symbol	SALL1
Entrez Gene	6299(Human)
SwissProt	Q9NSC2(Human)

\*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact [info@arexbio.com](mailto:info@arexbio.com) or your local distributor.

\*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

**DATASHEET****SALL1 Rabbit Polyclonal Antibody**

CAT. NO. APA18239

**DATA**

Western blot analysis of SALL1 expression in Y79 (A) whole cell lysates. (Predicted band size: 140 kD; Observed band size: 130 kD)

**STORAGE**

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

**NOTE**

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.