

DATASHEET

CCDC112 Rabbit Polyclonal Antibody

CAT. NO. APA18277

KEY FEATURES

Target	CCDC112	Source / Host	Rabbit
Reactivity	Human	Clonality	Polyclonal
Applications	WB	Conjugation	Unconjugated
Form / Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.		Storage at -20°C

BACKGROUND

Coiled-coil domain-containing protein 112 (CCDC112), also known as Mutated in bladder cancer protein 1. Subcellular location: Cytoplasm, cytoskeleton, microtubule organizing center, centrosome, centriolar satellite

APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:500 - 1:1000
----	----------------

*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

PRODUCT OVERVIEW


Description	Rabbit polyclonal antibody to CCDC112
Specificity	Recognizes endogenous levels of CCDC112 protein.
Antibody Type	Primary antibody
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the N-terminal region of human CCDC112. The exact sequence is proprietary.
Purification	The antibody was purified by immunogen affinity chromatography.
Molecular Weight	Predicted: 53 kD; Observed: 55 kD
Form/Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Alternative Names	MBC1; Coiled-coil domain-containing protein 112; Mutated in bladder cancer protein 1
Gene Symbol	CCDC112
Entrez Gene	153733(Human)
SwissProt	Q8NEF3(Human)

*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact info@arexbio.com or your local distributor.

*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

DATASHEET**CCDC112 Rabbit Polyclonal Antibody**

CAT. NO. APA18277

DATA

Western blot analysis of CCDC112 expression in ZR751 (A) whole cell lysates. (Predicted band size: 53 kD; Observed band size: 55 kD)

STORAGE

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

NOTE

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.