

**DATASHEET**

**ZNF12 Rabbit Polyclonal Antibody**

CAT. NO. APA18652

**KEY FEATURES**

Target	ZNF12	Source / Host	Rabbit
Reactivity	Human	Clonality	Polyclonal
Applications	WB	Conjugation	Unconjugated
Form / Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.		Storage at -20°C

**BACKGROUND**

Transcriptional repressor which suppresses activation protein 1 (AP-1)- and serum response element (SRE)-mediated transcriptional activity.

**APPLICATION**

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:500 - 1:1000
----	----------------

\*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

**PRODUCT OVERVIEW**


Description	Rabbit polyclonal antibody to ZNF12
Specificity	Recognizes endogenous levels of ZNF12 protein.
Antibody Type	Primary antibody
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the N-terminal region of human ZNF12. The exact sequence is proprietary.
Purification	The antibody was purified by immunogen affinity chromatography.
Molecular Weight	Predicted: 81 kD; Observed: 81 kD
Form/Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Alternative Names	GIOT3; KOX3; ZNF325; Zinc finger protein 12; Gonadotropin-inducible ovary transcription repressor 3; GIOT-3; Zinc finger protein 325; Zinc finger protein KOX3
Gene Symbol	ZNF12
Entrez Gene	7559(Human)
SwissProt	P17014(Human)

\*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact [info@arexbio.com](mailto:info@arexbio.com) or your local distributor.

\*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

**DATASHEET****ZNF12 Rabbit Polyclonal Antibody**

CAT. NO. APA18652

**DATA**

Western blot analysis of ZNF12 expression in NCIH292 (A) whole cell lysates. (Predicted band size: 81 kD; Observed band size: 81 kD)

**STORAGE**

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

**NOTE**

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.