

**DATASHEET**

**HSD17B7 Rabbit Polyclonal Antibody**

CAT. NO. APA18910

**KEY FEATURES**

Target	HSD17B7	Source / Host	Rabbit
Reactivity	Human	Clonality	Polyclonal
Applications	WB, FC	Conjugation	Unconjugated
Form / Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.		Storage at-20°C

**BACKGROUND**

Bifunctional enzyme involved in steroid-hormone metabolism and cholesterol biosynthesis . Catalyzes the NADP(H)-dependent reduction of estrogens and androgens and regulates the biological potency of these steroids. Converts estrone (E1) to a more potent estrogen, 17beta-estradiol (E2) . Converts dihydrotestosterone (DHT) to its inactive form 5a-androstane-3b,17b-diol . Converts moderately progesterone to 3beta-hydroxypregn-4-ene-20-one, leading to its inactivation . Additionally, participates in the post-squalene cholesterol biosynthesis, as a 3-ketosteroid reductase .

**APPLICATION**

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:500 - 1:1000
FC	1:10 - 1:50

\*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

**PRODUCT OVERVIEW**

Description	Rabbit polyclonal antibody to HSD17B7
Specificity	Recognizes endogenous levels of HSD17B7 protein.
Antibody Type	Primary antibody
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the N-terminal region of human HSD17B7. The exact sequence is proprietary.
Purification	The antibody was purified by immunogen affinity chromatography.
Molecular Weight	Predicted: 38 kD; Observed: 45 kD
Form/Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Alternative Names	3-keto-steroid reductase; 17-beta-hydroxysteroid dehydrogenase 7; 17-beta-HSD 7; Estradiol 17-beta-dehydrogenase 7
Gene Symbol	HSD17B7
Entrez Gene	51478(Human)
SwissProt	P56937(Human)

\*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact [info@arexbio.com](mailto:info@arexbio.com) or your local distributor.


\*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

**DATASHEET**


**HSD17B7 Rabbit Polyclonal Antibody**

CAT. NO. APA18910

**DATA**



Western blot analysis of HSD17B7 expression in 293 (A), HL60 (B), Jurkat (C) whole cell lysates. (Predicted band size: 38 kD; Observed band size: 45 kD)



Flow cytometric analysis of HL60 cells using Anti-HSD17B7 Antibody. The cells were fixed with 2% paraformaldehyde (10 min) and then permeabilized with 90% methanol for 10 min. The cells were incubated in 2% bovine serum albumin to block non-specific protein-protein interactions followed by the antibody at 37 °C for 60 min. The secondary antibody Goat Anti-Rabbit IgG (H&L) - AREX® Fluor 488 was incubated at 37 °C for 40 min. Isotype control antibody (blue line) was used under the same condition.

**STORAGE**

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

**NOTE**

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.